

Comprehensive rescue system



Technical Specifications

Empowering Bystanders to Rescue with Confidence

When medical emergencies occur, the faster care can be provided, the better the outcome. Mobilize Rescue Systems are the only trauma and medical systems that include the medical supplies and interactive instructions bystanders need to respond to medical emergencies in the critical minutes before EMS arrives.

Mobilize Rescue Systems are designed to help untrained bystanders assess, manage, and monitor a spectrum of medical emergencies according to the protocols and best practices of the American Heart Association (AHA), Tactical Emergency Casualty Care (TECC), and the American College of Surgeons.

Portable Trauma and Emergency Medical Response Kit

More than a first aid kit or a stop-the-bleed kit, the Comprehensive rescue system includes medical supplies to treat a multitude of life-threatening injuries and medical issues such as severe bleeding, airway obstructions, traumatic injuries, allergic reactions, and broken bones. It also includes space for additional non-included supplies such as EpiPens and Narcan.*

Integrated App Provides Real-time Instructions

The Comprehensive rescue system features the integrated Mobilize Rescue application, accessible through a built-in tablet in the unit. Step-by-step instructions on how to use the supplies are provided through simple interactive prompts, written instructions, and visual cues that guide rescuers to provide care during a medical emergency.

Peace of Mind that You Are Prepared

Preparing for cardiac emergencies, traumatic injuries, and life-threatening medical issues is key for safety leaders across all industries. By installing the Comprehensive rescue system, safety leaders can be confident, knowing they are prepared – and are preparing others – to help during critical medical emergencies that can become fatal during the average 9-15 minutes it takes for an ambulance to arrive.¹

* Prescription medications are not included in Mobilize Rescue Systems.

System Contents

Trauma Supplies

SOF-T Wide Tourniquet	4
QuikClot Bleeding Control Dressing	2
6" NAR Emergency Trauma Dressing	2
Hyfin Chest Seal	2
Water-Jel Universal Burn Dressing	2
Triangular Bandage	2
4.5" Sterile Conforming Stretch Gauze	2
5" x 9" Sterile Combine ABD Pads	4
10" x 30" Sterile Multi-Trauma Dressing	2
36" SAM Emergency Splint	1
4" Elastic Wrap Bandage	2
4" x 5" Cold Compress	1
Adhesive Tape (2.5 yd.)	1

Medical Supplies

CPR Face Shield with Bite Block	1
81mg Chewable Aspirin (Bottle)	1
12mg Dissolvable Allergy Tablets (Box)	1
15mg Insta-Glucose	1
Emergency Space Mylar Blanket	2

Other

Portable Charger* and Charging Cord	1
USB Charging Cube	1
Inspection Card	1
Proof Seals	12
Bag with Biohazard Markings	2
Trauma Shears	1
Nitrile Gloves	10

OSHA-Compliant First Aid Supplies

Eye Pads	2
Adhesive Bandages (Assorted Sizes)	16
Burn Cream Packets	10
Triple Antibiotic Ointment Packets	10
Antiseptic Wipes	2
Tweezers	1
Hand Sanitizer Packets	10
Eye Wash	1
First Aid Guide	1

Rescue Application

Multi-user app license

*5 Volt Charger – for temporary backup power in the event of battery depletion on the tablet.

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824, USA
978-421-9655
800-804-4356

For subsidiary addresses and fax numbers as well as other global locations, please go to www.zoll.com/contacts.

Comprehensive Rescue System

Equipment


Portable trauma and emergency medical response kit with a built-in Android tablet and interactive instructional rescue app.

The system's durable, water-resistant case provides maximum protection for supplies in a range of environments and conditions.

Dimensions: | 17.2" X w 13.7" X d 6.5"

Weight: 15 lbs

Backpack Version Available



10" Application Display – Android LCD tablet

Technology

The Mobilize Rescue app is accessed through the tablet within the unit. The app provides real-time treatment instructions through interactive text and visual prompts.

App Integration

All the medical supplies are clearly labeled, logically organized, and color-coded to match the instructions in the Mobilize Rescue app. This system's intuitive organization makes locating supplies easy, saving valuable time during an emergency.



Alarmed wall cabinets available.

Treatment Capabilities

SEVERE BLEEDING	OPIATE OVERDOSES	AMPUTATIONS	TRAPPED LIMBS
UNCONSCIOUSNESS	HYPOTHERMIA	ALLERGIC REACTIONS	BURNS
CPR INSTRUCTIONS	SEIZURES	CONFUSION	FRACTURES & SPRAINS
AED APPLICATION	CHOKING	STINGS	DIABETIC PROBLEMS
TROUBLE BREATHING	IMPALED OBJECTS	PAIN	CHEST PAIN

STORAGE RECOMMENDATIONS & DIAGNOSTIC CHECK | Individual supplies within the unit may be negatively impacted by prolonged exposure to hot and cold temperatures. Individual supplies have expiration dates that vary from 1-5 years from date of purchase; it is crucial that you perform diagnostic checks to ensure that the unit is charged to 100%, the proof seal is intact (not tampered with), and that the medical supplies are not damaged or contaminated.

The Comprehensive Unit includes a tablet device that should be used only in the temperature range of 0°C (32°F) to 40°C (104°F) (storage -20°C (-4°F) to 60°C (140°F)) to avoid damage.

¹Journal of Emergency Medical Services. (2005.) Response Times: Myths, Measurement and Management. Retrieved from: <https://www.jems.com/articles/2005/08/response-times-myths44-measure.html>
Specifications subject to change without notice.

©2020 ZOLL Medical Corporation. All rights reserved. Mobilize Rescue Systems and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

MCN PP 2001 0404

